




A SUBSIDIARY OF VITAE CORE COMPANY



An ISO 9001-2015 Company

VGRIPP  
**NOVA**  
Variceal Ligation System



Used for treatment of esophageal varices above the GE junction the system incorporates automatic alerts for band release and count of unused bands.



6 Latex rubber bands in black colour

helps easy visualization

Excellent 100% vision

ensure total viewing

Minimum release force required for band

causing least compression of the endoscope

In situ count of unreleased bands

facilitates proper planing

Special handle adapters

keeps the handle secure and firm

Convenient irrigation port for flushing

blood or any fluids which obstruct vision

Single suture release

keeps the lumen of the channel free for the varix and vision

Catheter loaded trip cord

ensures insertion of trip cord in one pass

## Product Specifications

Ref No.	111110
Item	VGRIPP NOVA Variceal Ligation System
Contents	Ligation cylinder with 6 bands, release handle, release cord with channel cover NOVA, 2 nos release handle adaptors
Ref No.	111120
Item	VGRIPP NOVA Variceal Ligation Cylinder
Contents	Ligation cylinder with 6 bands, release cord with channel cover NOVA
Ref No.	111130
Item	VGRIPP NOVA Release handle
Contents	Release handle, 2 nos release handle adaptors

Regd. Office - 344/5, N.S.C. Bose Road, Kolkata - 700 047, India. Tel - (033)2421 5614

Website : [www.imi.co.in](http://www.imi.co.in)



A SUBSIDIARY OF VITAE CORE COMPANY



# VGRIPP

Vgripp Variceal Ligation Cylinder



Used for treatment of esophageal varices above the GE junction the system incorporates automatic alerts for band release and count of unused bands.



# VGRIPP

Vgripp Variceal Ligation Cylinder

6 Latex rubber bands in black color

helps easy visualization

Excellent 100% vision

ensures total viewing

Minimum release force required

causing least compression of the endoscope

In situ count of unreleased bands

facilitates proper planning

Audible indicator

for the release of varices and target the next varix

Single suture release

keeps the lumen of the channel free for the varix and vision

Catheter loaded trip cord

ensures insertion of trip cord in one pass

## Product Specifications

Ref No.	111210
Item	VGRIPP Variceal Ligation Cylinder
Contents	Ligation cylinder with 6 bands, release cord with channel cover VG

Ref No.	111220
Item	VGRIPP Release Handle
Contents	Release handle with elastic band

Regd. Office - 344/5, N.S.C. Bose Road, Kolkata - 700 047, India. Tel - (033)2421 5614

Website : [www.imi.co.in](http://www.imi.co.in)



A SUBSIDIARY OF VITAE CORE COMPANY



An ISO 9001-2015 Company

# VGRIPP

Multiband Variceal Ligation Cylinder

6 Latex rubber bands in black colour  
Double suture release band  
Loading catheter  
A stopper to ensure better suction



Used for the treatment of esophageal varices above the GE junction

## Product Specifications

Ref No.	111500
Item	VGRIPP Multiband Variceal Ligation Cylinder
Contents	Ligation cylinder with 6 bands and release cord (1), loading catheter

Recommended Scope Diameter – 9.5 mm – 11 mm

Regd. Office - 344/5, N.S.C. Bose Road, Kolkata - 700 047, India. Tel - (033)2421 5614

Website : [www.imi.co.in](http://www.imi.co.in)



A SUBSIDIARY OF VITAE CORE COMPANY



AN ISO 9001-2015 Company

# PEG KIT

Percutaneous Endoscopic  
Gastroscopy



The PEG KIT is placed percutaneously for nutrition and medications in patients with dysphagia or neurological disorders.

The peg tube helps in the transport of food and other liquids

The feeding adapter helps for the inlet of medicine and other fluids into the peg tube


The bolster is used to affix the tube firmly

The blue pull wire, scalpel and introducer needle helps in the procedural set up of peg tube

## Product Specifications

Ref. No. 701110  
Item EzFeed Peg Set 20F (Full Kit)  
Contents BLACKWIRE PULL WIRE, SEMI KIT COMP EXCLUDING, BLUEWIRE, SNARE, HEMOSTATS, SCISSORS, GAUZE 5PCS, DRAPE 1PC, SYRINGE 5ML-1PC, NEEDLE 22G & 25G, BOLUS, TWIST LOCK.

Ref. No. 702110  
Item EzFeed Peg Set 24F (Full Kit)  
Contents BLACKWIRE PULL WIRE, SEMI KIT COMP EXCLUDING, BLUEWIRE, SNARE, HEMOSTATS, SCISSORS, GAUZE 5PCS, DRAPE 1PC, SYRINGE 5ML-1PC, NEEDLE 22G & 25G, BOLUS, TWIST LOCK.



# PEG

## Percutaneous Endoscopic Gastroscopy

The PEG KIT is placed percutaneously for nutrition and medications in patients with dysphagia or neurological disorders.



### Product Specifications

Ref. No. 702100  
Item EzFeed Peg Set 24 Fr (Mini Kit)  
Contents 1 No. each of pack tube, feeding adapter, bolster, scalpel, clamp, introducer needle and blue pull wire.

Ref. No. 701100  
Item EzFeed Peg Set 20 Fr (Mini Kit)  
Contents 1 No. each of pack tube, feeding adapter, bolster, scalpel, clamp, introducer needle and blue pull wire.

**EzFeed Peg Replacement Gastrostomy Tube**  
used as replacement gastrostomy tube through an established gastrostomy tract.

Ref. No. 712400  
Item EzFeed Peg Replacement Tube 24 Fr

Ref. No. 713400  
Item EzFeed Peg Replacement Tube 26 Fr



A SUBSIDIARY OF VITAE CORE COMPANY



An ISO 9001-2015 Company

A SUBSIDIARY OF VITAE CORE COMPANY



An ISO 9001-2015 Company



**Regd. Office** - 344/5 N.S.C. Bose Road. Kolkata - 700 047, India

Ph. No. : (033)2421 5614

**E-mail** : [info@imi.co.in](mailto:info@imi.co.in)

**Website** : [www.imi.co.in](http://www.imi.co.in)